STERILE R

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INGLÉS

DESCRIPTION AND ACTION FOUNDATIONS

Bonefill comprises a mineral structure removed from the bone of bovine femur, that presents a structure similar to human bone in its inorganic portion, free of organic contaminations and it is a product with high conductive potential, an implantable medical/dental biomaterial, which is used in the reconstruction of bone losses where the bone neoformation or bone remodeling is expected.

It is produced from the natural bone extracted from bovine femur, completely deproteinised. The fresh bone is crushed and submitted to a sequence of baths that solubilise the organic structures such as cells remaining from organic matrix, fibers and proteins, remaining only the mineral portion and avoiding the induction of possible immunogenic processes to the organism.

Bonefill is used as material for bone grafting, for the replacement or neoformation of bone tissue, periodontal bone defects, points of teeth extraction and other bone defects at the mouth, maxillofacial, base of skull, spinal fusion areas, besides orthopedic areas. It can be also used as bone graft around implants and as a support in the treatment of maxillary sinus elevation. The expected time for incorporating the products made of mineralized bovine bone is of 120 to 180 days. It is a biocompatible material presented in proper granulation and with bone-conducting characteristics helping the guided tissue regeneration processes, periodontal defects and bone neoformation. The granulated Bonefill is delivered in glass vial containing a deproteinised powdered product, whilst the blocks are presented in blister and tyveck.

PRODUCT'S COMPOSITION

Bonefill is manufactured from natural bone extracted from bovine femur, totally deproteinised

INDICATIONS AND PURPOSE OF USE

Bonefill – material for bone grafting, which is used as material for bone grafting for replacement or neoformation of bone tissue, periodontal bone defects, points of teeth extraction and other bone defects at the mouth, maxillofacial, base of skull, spinal fusion areas, besides orthopedic areas. It can be also used as bone graft around implants and as a support in the treatment of maxillary sinus elevation and bone neoformation.

The incorporation of the bone graft is accomplished by neoformation of bone that comes from the recipient bed within the period from 120 to 180 days.

It is supplied in the sterile condition, since it is kept under ideal conditions of storage and conservation and that the integrity of the package is not jeopardized. It is sterilized through

Gamma Radiation (25 kGy) and it should not be used in case the validity term is expired. SOLE USE – Do not re-sterilize or reprocess and destroy after the use.

PRECAUTIONS, RESTRICTIONS, WARNINGS

1. STERILE - as long as maintained the integrity of the packaging, period of validity and storage conditions.

2. Professional use only is the responsibility of the dentist or doctor their prior training to use this product. Only qualified professionals with expertise in surgical techniques and procedures necessary for proper use of the product should make use of Bonefill. The use of incorrect surgical techniques may cause discomfort as painful sensation, hypoesthesia, edema. The size of the particles should not be changed.

3. PROHIBITED RESTERILIZE AND REPROCESS - If re-sterilized or reprocessed may occur change at the physical chemical properties and crystallinity levels of Bonefill causing foreign body reaction.

4. PROHIBITED REUSE - If reused or used with expired validity, may cause irritation, infection, inflammation and other adverse events, compromising the health and safety of the patient. Bionnovation does not recommend reuse, re-processing or re-sterilization, discard it as current legislation for medical waste.

5. The use of the product with surgical techniques and inadequate biosecurity conditions may damage the patient leading to unsatisfactory results

6. Always sterilize the tools before using them.

7. The clinical and radiographic evaluation must be done prior to surgery, to help the correct treatment planning. Determination of bone quality and quantity, repairs and anatomical structures and analysis of neighboring teeth.

8. In all surgeries involving Bonefill particles must be observed proper techniques used for asepsis and antisepsis.

9. The abuse of alcohol, tobacco, drugs, steroids or lack of proper oral hygiene can significantly impair the success of the treatment.

10. All potential adverse effects as dehiscence, inflammation, bone loss, hemorrhage, must be previously informed to the patient.

11. It is supplied in sterile condition and once opened should be used in aseptic conditions. Should always work with sterile fields, appropriate instruments to the procedure and in good condition in order to eliminate sources of infection and damage to the product.

12. Bonefill should be used only for the purpose for which it is intended.

13. If occur complications impossible to be controlled, tissue inflammation or evidence of infection is recommended the immediate removal of the material.

14. Bonefill is provided in sterile double packaging (25 kGy gamma radiation). Provided that the package integrity is not compromised in any way, it will save the sterile product up to 3 years from the date of sterilization.

15. There are no restrictions as to maximum amount of product that can be deployed. The amount will be determined by the professional after analyzing the size of the surgical site. 16. The surgeon should evaluate the indication in patients who are carriers of diseases or are making use of medication that may alter the repair metabolism.

17. The remainder of the packaging material should not be reused, reprocessed or re-sterilized, discard it mischaracterized as current legislation for medical waste, do not discard contaminated products in household waste.

18. In cases of adverse effects occurring in patients, the responsible must contact immediately with the SAC Bionnovation (Customer Service) by the number 0800 707 3824 or by e-mail sac@bionnovation.com.br. The Bionnovation Biomedical Products is responsible for notifying the ANVISA (Health Surveillance Agency) about the relevant occurrences according to internal technovigilance procedure.

19. Bonefill was developed in order to prevent that its use does not compromise the clinical condition of patients as well as their safety.

20. It is recommendable that the surgeon keeps a follow-up report on the progress of the case together with the results.

21. All Bonefill bone replacements of bovine origin are manufactured with bones of animals from the Brazilian herd tracked by the SISBOV system. As per the geographic risk regulations issued by the International Zoonosary Code and the Scientific Steering Committee of the European Union (SSCEC of August 2005), Brazil is Bovine Spongiform Encephalopathy (BSE)-free. However, according to Administrative Rule 516/97 even Brazil declaring itself Bovine Spongiform Encephalopathy-free, the processing to which the products are submitted are known to effectively inactivate the BSE agent, and the animals used to manufacture Bonefill products are registered in the Brazilian bovine and bubaline traceability system - SISBOV, every product of bovine origin, even remotely, carries a risk of transmitting BSE. In the event Creutzfeldt - Jakob disease (CJD) symptoms appear, the professional must notify the health authorities.

22. The product Bonefill - Materials for Bone Graft is biocompatible, has no cytotoxicity, acute systemic toxicity, subchronic toxicity, carcinogenicity, genotoxicity, and is not a sensitizing product.

We recommend that the identification stickers that come with the product be attached to the patient's documentation: clinical record of the patient, the report delivered to the patient, product sales invoice, vendor control and control of the surgeon in charge,

ensuring the complete traceability of the product.

Contra indications

1. Bonefill, as well as all the other biomaterials, should not be placed on an existing active infection or in any other degenerative disease that might affect the biomaterial's placement.

2. It must not be utilized in patients that are notable, under the clinical point of view, to be submitted to a medical or odontological intervention. Such as, for example, in patients with uncompensated diabetes.

3. Bonefill is not indicated for odontopediatric patients.

4. It's contraindicated for procedures different from those recommended in item "Use Indication".

5. Bonefill should not be exposed to the external medium.

STERILIZATION

Bonefill is supplied STERILE (Gamma Radiation 25kGy), provided the package integrity is not compromised.

PRE AND POST-SURGICAL CARES

The pre-surgical evaluation, the correct indication of materials and the employment of compatible techniques and procedures, as well as the post-surgical follow-up and controls, are fundamental to achieve the desired results.

Pre-Surgical: All the patients that will be submitted to a surgical procedure must be carefully examined and evaluated with the purpose of determining their clinical and radiographic state, as well as their dental, bone, or adjacent soft tissue deficits that might influence the final result of the intervention.

Post-Surgical: Please, note the post-surgical cares for the surgical procedures. Painkillers, antibiotics, or rest for 24-48 hours may be prescribed, varying as a function of the procedure and of the professional technical conduct. The product must not be exposed to the mouth environment after the immediate post-operative. There must be a good coaptation of the surgical snip borders, in order not to be any contamination, which would compromise the surgeon's result. Exposure to the mouth environment drastically reduces the absorption time.

SPECIAL CONDITIONS FOR THE PRODUCT'S STORAGE AND TRANSPORTATION, CONSERVATION AND/OR HANDLING

Storage and Transportation:

Store and transport the product away from direct sunlight, and from heat. (temperature between 15-25°C) and humidity sources. Keep the packaging sealed until its utilization time. Do not use it if the sterile package has been opened, or if it's damaged, or if the sterilization validity date has expired in order to avoid possible contamination. Discard any mischaracterized product according to the applicable legislation for hospital waste, or return the damaged packages to the factory, including the device.

Conservation and Handling

In case of any alteration in the Bonefill characteristic that might have mischaracterized it, please discard it according to the applicable legislation for hospital waste, or return the damaged packages to the factory, including the device.

COMMERCIAL PRESENTATION FORMS

Bonefill available in the form of granules and blocks.

Bonefill with fine granulation (particles with a diameter of up to 0.6 mm), average granulation (particles with diameter between 0.6 mm and 1.5 mm), and large granulation (particles with diameter between 1.5 mm and 2.5 mm)

USE INSTRUCTIONS

Bonefill granules

1. Observe the preparation of the site and use of proper sterile garments by the surgeon, assistants and patient for a sterile procedure.

2. Before opening the package, a vessel for hydration of Bonefill should be prepared and association with other grafting materials (autogenous bone tissue collected during the surgery so as to improve the bone neoformation process). The vessel should be sterile, as well as the physiological saline that is used for hydration.

3. Remove the glass vial from the final and secondary package.

4. Using sterile gloves open the primary blister only at the site prepared for the surgical procedure.

5. Remove the seal and the lid.

6. Add the vial content in a sterile vat.

7. Mix sterile saline solution with the bone graft with a sterile spatula, the product will be

increased in volume due to hydration.

8. Before applying the product, the bone area that will receive the biomaterial should be opened and scraped with a curette, decorticated or perforated for exposure of organic matrix and the whole jeopardized tissue should be removed. Conventional instruments such as currettes, spatulas, plastic or metallic applicators can be used so as to ease the application at the surgical site.

9. Irrigate the site before putting the graft

10. Make sure there is proper hemostasis to ensure that the material is retained in the surgical site.

11. Use the product immediately after reconstitution with physiological solution.

12. Ensure the maximum contact between the bone replacement material and the receiving bone, with the total filling of space and good compaction.

13. Reposition the piece on the grafted area and ensure that the coverage of the surgical bed is complete. There cannot be the exposure of biomaterial, because it is basic for the surgical success and to avoid contamination of the area.

Note: The professional shall evaluate the characteristics of the proper bone support for introduction, once the bone density is a determinant and a variable factor when planning and conducting endosseous implants. Every type of density presents several advantages and disadvantages, which will determine the most proper planning. The bone density impacts on surgical technique and repair time. If such is the case, the previous knowledge on bone density becomes an essential factor for a correct planning.

Bonefill as block

1. Observe the preparation of the site and use of proper sterile garments by the surgeon, assistants and patient for a sterile procedure.

2. Before opening the package, a vessel for hydration of Bonefill should be prepared and association with other grafting materials (autogenous bone tissue collected during the surgery so as to improve the bone neoformation process). The vessel should be sterile, as well as the physiological saline that is used for hydration.

3. Using sterile gloves open the primary blister only at the site prepared for the surgical procedure.

4. Add the vial content in a sterile vat.

5. Mix sterile saline solution with the bone graft with the aid of sterile tweezers, and the product will be increased in volume due to hydration.

6. Cut out the block with the aid of a sterile disc according to the needs of the patient.

7. Before applying the product, the bone area that will receive the biomaterial should be

opened and scraped with a curette, decorticated or perforated for exposure of organic matrix and the whole jeopardized tissue should be removed.

8. Irrigate the site before putting the graft.

9. Make sure there is proper hemostasis to ensure that the material is retained in the surgical site.

10. After the exposure of medullary portion of the bone, Bonefill as block can be attached with the aid of Bolt for graft and attachment, providing higher stabilization for bone block.

11. Ensure the maximum contact between the bone replacement material and the receiving bone, with the total filling of space.

12. Reposition the piece on the grafted area and ensure that the coverage of the surgical bed is complete. There cannot be the exposure of biomaterial, because it is basic for the surgical success and to avoid contamination of the area.

CARE WHEN DISCARDING THE PRODUCT

The product's disposal must comply with the environmental and biosafety laws in force. Do not discard contaminated products in the general waste.