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## DESCRIPTION AND ACTION BASIS

Surgitime PTFE is a synthetic non-absorbable membrane that is 100% bio-compatible and not derived from animal source. Considered as a barrier for tissue regeneration, it is indicated for regeneration procedures in Periodontology, Implant Dentistry, orthopedics, neurosurgery and other dentistry or medical surgeries, especially for reconstruction of tissues. Polytetrafluoroethylene (PTFE) membranes or mechanical barriers for Guided Tissue Regeneration (GTR) are used to prevent migration of cells from epithelial and connective tissues, what would cause bone growth inhibition, thus providing a proper space for the formation of a natural fibrin structure, which is the bone tissue precursor. The membrane provides a space between the flap and the bone tissue and its tissue isolating property promotes tissue growth.

## PRODUCT COMPOSITION

Polytetrafluoroethylene (PTFE)

## INDICATIONS AND USE PURPOSE

Surgitime PTFE is used in regenerative techniques of periodontics, implantology, orthopedics, neurosurgery, other dentistry or medical surgeries or any surgical procedure requiring a mechanical barrier, such as the treatment of horizontal and vertical periodontal defects, formation of new bone in alveolar ridges, protection against epithelial invagination in sinus lift procedures and formation of proximal areas around dental implants. It is available in different shapes for better adaptation to the receiving site and it can also be molded by using scissors or a sterile scalpel.

## PRECAUTIONS, RESTRICTIONS AND WARNINGS

- Do not use if the package is damaged or after the expiration date.
- Abuse of alcohol, tobacco, drugs and corticosteroids or inappropriate oral hygiene may significantly compromise the success of the treatment.
- Indistinguishably dispose or return the damaged or expired products to the manufacturing site.
- Do not re-sterilize, re-process or reuse. .
- The membrane should be used for the indicated purposes only.
- We recommend attaching the identification adhesive labels accompanying the product to the patient's records.
- The rest of the packing material should be disposed of as per the legislation in force for hospital waste; do not dispose contaminated products in ordinary waste. Destroy after use.
- Surgitime PTFE is supplied in sterile packages.
- As any other membrane, Surgitime PTFE should not be used in the event of preexisting active infection.

- In the event of complications that cannot be controlled, tissue inflammation or evidence of infection, the immediate removal of the material is recommended.
- In all procedures involving the use of Surgitime PTFE, appropriate asepsis and antiseptic techniques should be observed.
- Do not use in pediatric patients.
- Surgitime PTFE is supplied in sterile packages (ETO Radiation). Providing the packaging integrity is preserved, it will keep the product sterile up to 2 years from the sterilization date.
- Surgitime PTFE should not be used in patients that are not eligible, from the clinical point of view, to undergo an oral surgery.
- There is no restriction regarding the maximum amount of product that can be implanted. The amount should be determined by the professional performing the procedure after analysis of the surgical bed.
- The packaging should be opened just before use.
- It is supplied in sterile form and after the package is opened, it should be used under aseptic conditions.
- Sterile fields, suitable instruments for the procedure and in good conditions of use should be used in order to eliminate infection sources and damages caused to the components by inappropriate instruments.
- Only qualified professionals experienced in surgical techniques and in the required procedures for the proper use of this product should use Surgitime PTFE.
- The product should be used under aseptic conditions and employing sterile compatible instruments.
- Do not use if the package is damaged or after the expiration date.
- Store in a clean place away from heat and humidity sources.
- Clinical and radiological evaluation should be made before the implantation surgery in order to help in the correct planning of treatment. Determination of bone amount and quality, repairs and anatomical structures and analysis of adjacent teeth is also important.
- In the event of adverse events in patients, the professional in charge should contact Bionnovation Consumer Services (SCA) immediately at 0800 707 3824 or e-mail sac@bionnovation.com.br. Bionnovation Produtos Médicos is responsible for notifying ANVISA (National Agency for Sanitary Surveillance) about relevant occurrences and the pertinent procedures are described in the internal rules for Quality General Procedures – "Procedures to Notify Health Authorities about Incidents".

STERILE.  
FOR PROFESSIONAL USE ONLY.  
NON-REUSABLE PRODUCT.

ALWAYS STERILIZE THE INSTRUMENTS BEFORE USING THEM.

## STERILITY

Surgitime PTFE is supplied as a STERILE product (Ethylene Oxide - ETO). Providing the package integrity is kept.

## POST-SURGERY CARE

Painkillers, antibiotics and rest during the first 24-48 hours can be prescribed after surgery, varying according to the procedure and the professional technical conduction. The product should not be exposed to the oral medium immediately after surgery. The surgical flap edges should be properly co-adapted in order to prevent exposition, what would compromise the surgery result. Exposure to the oral medium may cause the accumulation of bacterial plaque on the membrane surface.

## STORAGE AND TRANSPORTATION

The product should be transported and stored away from direct sunlight, heating sources or humidity. Keep the package sealed until its use. Discard the product if package is damaged or after the expiration date.

## PRODUCT PRESENTATION

Surgitime PTFE is a non-absorbable barrier made with polytetrafluoroethylene, with thickness of 0,1mm, 0,25mm and 0,50mm. It is supplied in different shapes: Surgitime PTFE "H", Surgitime PTFE Trapeze, Surgitime PTFE 3,0 x 2,0 cm Rectangle, Surgitime PTFE 6,0 x 6,0 cm, Surgitime PTFE 8,0 x 8,0 cm, Surgitime PTFE 10,0 x 4,0 cm Rectangle, Surgitime PTFE 10,0 x 6,0 cm Rectangle, Surgitime PTFE 10,0 x 8,0 cm Rectangle, Surgitime PTFE 10,0 x 10,0 cm, Surgitime PTFE 12,0 x 6,0 cm Rectangle, Surgitime PTFE 12,0 x 12,0 cm.

## INSTRUCTIONS

1. Place the package contents on a sterile surgical field.
2. Detach the flap so that later the membrane exceeds at least 2 mm the area to be protected.
3. Using the aseptic surgery techniques applicable in the case, prepare the receiving bed for the membrane.
4. Cut the membrane in a suitable size, with the aid of a sterile pair of scissors, aiming maximum adaptation.
5. Adapt the membrane to the site, observing its margins. All of it should be over soft tissue and without folds.

6. Re-place the flap over the membrane.
7. Suture without involving the membrane. Only in cases of use for medical field, the membranes with thickness of 0,1mm, 0,25mm and 0,50mm can be involved.
8. Using surgical cement is facultative to the surgeon in charge.
9. Antibiotics, painkillers or anti-inflammatory drugs may be used after surgery.
10. The membrane can be removed once its purpose is accomplished.